



dualX transforms the fusion environment by delivering a powerful implant through a minimally-invasive approach

Please send your completed scrub sheets to:
customerservice@amplifysurgical.com
 765-AMPLIFY (765-267-5439)

Top View



Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Lateral Lumbar Interbody Fusion (LLIF)								
1283-1005-002	dualX LLIF 7-9mm x 45mm	0°	13	22	7	9	54	45
1283-1005-003	dualX LLIF 7-9mm x 50mm	0°	13	22	7	9	59	50
1283-1005-004	dualX LLIF 7-9mm x 55mm	0°	13	22	7	9	64	55
1283-1007-001	dualX LLIF 8-11mm x 40mm	7°	13	22	8	11	49	40
1283-1007-002	dualX LLIF 8-11mm x 45mm	7°	13	22	8	11	54	45
1283-1007-003	dualX LLIF 8-11mm x 50mm	7°	13	22	8	11	59	50
1283-1007-004	dualX LLIF 8-11mm x 55mm	7°	13	22	8	11	64	55
1283-1007-012	dualX LLIF 10-13mm x 45mm	7°	13	22	10	13	54	45
1283-1007-013	dualX LLIF 10-13mm x 50mm	7°	13	22	10	13	59	50
1283-1007-014	dualX LLIF 10-13mm x 55mm	7°	13	22	10	13	64	55
1283-1007-022	dualX LLIF 12-15mm x 45mm	7°	13	22	12	15	54	45
1283-1007-023	dualX LLIF 12-15mm x 50mm	7°	13	22	12	15	59	50
1283-1007-024	dualX LLIF 12-15mm x 55mm	7°	13	22	12	15	64	55
1283-1000-006	dualX LLIF 10-13mm x 40mm	12°	13	22	10	13	49	40
1283-1000-008	dualX LLIF 10-13mm x 50mm	12°	13	22	10	13	59	50
1283-1000-009	dualX LLIF 10-13mm x 55mm	12°	13	22	10	13	64	55
1283-0009-001	dualX LLIF Lock Screw 2 - 40mm				N/A			
1283-0009-002	dualX LLIF Lock Screw 2 - 45mm				N/A			
1283-0009-003	dualX LLIF Lock Screw 2 - 50mm				N/A			
1283-0009-004	dualX LLIF Lock Screw 2 - 55mm				N/A			

Note: dualX expansion in width and height is definitive and not variable. For example, the dualX LLIF 0 DEG x 7-9mm x 45mm expands exactly from 13mm to 22mm (9mm expansion) in width and 7mm to 9mm (2 mm) in height.

Please refer to package insert for current warnings, precautions and instructions for use.
 Copyright 2021. All Rights Reserved. dualX is a registered trademark of Amplify Surgical, Inc.
 For more information call 765-AMPLIFY, email info@amplifysurgical.com or visit amplifysurgical.com
 ML10014.E



Top View




Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal Lumbar Interbody Fusion (TLIF)								
1253-1006-001	dualX TLIF 7-9mm	0°	12	21	7	9	36	30
1253-1006-004	dualX TLIF 8-11mm	8°	12	21	8	11	36	30
1253-1006-005	dualX TLIF 9-12mm	8°	12	21	9	12	36	30
1253-1006-006	dualX TLIF 10-13mm	8°	12	21	10	13	36	30
1253-1006-007	dualX TLIF 11-14mm	8°	12	21	11	14	36	30
1253-1006-008	dualX TLIF 12-15mm	8°	12	21	12	15	36	30
1253-0033	dualX TLIF Lock Screw 2	N/A						

Top View



Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal/Posterior Lumbar Interbody Fusion (T/PLIF)								
1282-1000-001	dualX T/PLIF 7-9mm	0°	12	21	7	9	36	22
1282-1001-014	dualX T/PLIF 8-10mm	8°	12	21	8	10	36	22
1282-1000-004	dualX T/PLIF 8-11mm	8°	12	21	8	11	36	22
1282-1000-005	dualX T/PLIF 9-12mm	8°	12	21	9	12	36	22
1282-1000-006	dualX T/PLIF 10-13mm	8°	12	21	10	13	36	22
1282-1000-007	dualX T/PLIF 11-14mm	8°	12	21	11	14	36	22
1282-1001-016	dualX T/PLIF 10-12mm	12°	12	21	10	12	36	22
1282-1001-017	dualX T/PLIF 10-13mm	12°	12	21	10	13	36	22
1282-1001-018	dualX T/PLIF 11-13mm	15°	12	21	11	13	36	22
1282-1000-019	dualX T/PLIF 11-14mm	15°	12	21	11	14	36	22
1282-0005	dualX T/PLIF Lock Screw 2	N/A						