

Please send your completed scrub sheets to:
customerservice@amplifysurgical.com
 765-AMPLIFY (765-267-5439)

dualXSlim transforms the fusion environment by delivering a powerful implant through a minimally-invasive approach

Top View



Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal/Posterior Lumbar Interbody Fusion (T/PLIF)								
1294-1001-004	dualXSlim T/PLIF 7-9mm	8°	10	20	7	9	35	23
1294-1001-005	dualXSlim T/PLIF 7-10mm	12°	10	20	7	10	35	23
1294-1001-012	dualXSlim T/PLIF 8-11mm	15°**	10	20	8	11	35	23
1294-1001-013	dualXSlim T/PLIF 9-12mm	15°**	10	20	9	12	35	23
1294-1001-022	dualXSlim T/PLIF 10-13mm*	15°**	10	20	10	13	35	23
1294-1001-030	dualXSlim T/PLIF 11-14mm*	15°**	10	20	11	14	35	23
1284-0006-001	dualXSlim Lockout Screw 2				N/A			

* Available upon request

** 18°, 20° Lordotic options available upon request

Note: dualXSlim expansion in width and height is definitive and not variable. For example, the dualXSlim 7-9mm x 8 DEG expands exactly from 10mm to 20mm (10mm expansion) in width and 7mm to 9mm (2 mm) in height.

Please refer to package insert for current warnings, precautions and instructions for use.

Copyright 2025. All Rights Reserved. dualX is a registered trademark of Amplify Surgical, Inc.

For more information call 765-AMPLIFY, email info@amplifysurgical.com or visit amplifysurgical.com

ML10014.F