


AMPLIFY<sup>®</sup>  
SURGICAL

# EXPAND THE POSSIBILITY



dualX<sup>®</sup>Slim

**Dual Expanding** Interbody Fusion System

# dualX<sup>®</sup>Slim

## Dual Expanding Interbody Fusion System

dualX<sup>®</sup>Slim transforms the fusion environment from insertion to spinal restoration by delivering a powerful dual-expanding implant through a minimally-invasive approach.




One Step Expansion

## Dual Expanding Interbody Implant Engineered for Spinal Fusion



Collapsed width designed to reduce neural retraction  
Low profile protects endplates during insertion

**10 mm Insertion Width**




Minimize subsidence with the largest footprint in the expandable cage market

Lateral expansion establishes stable footprint

Large, center bone graft chamber for post-expansion grafting

**20 mm Fully Expanded Width**



Powerful vertical expansion restores disc height for decompression


Vertical expansion assists in direct and indirect decompression

Multiple lordotic angles restore sagittal balance

**8°, 12°, 15°, 18°\***

**\*Available Upon Request**

Innovative dual locking mechanism designed to maintain integrity of implant until patient is fused



\*Available upon request



## Three Unique Implants for Posterior, Transforaminal and Lateral Approaches



### T/PLIF

- Heights 7mm\* expanding to 17mm\*
- Width 12mm expanding to 21mm
- 0°, 8°, 12° and 15° Lordosis\*



### TLIF

- Heights 7mm\* expanding to 17mm\*
- Width 12mm expanding to 21mm
- 0°, 8°, 12° and 15° Lordosis\*



### LLIF

- Heights 7mm\* expanding to 17mm\*
- Width 13mm expanding to 22mm
- Final Length 40 to 60mm
- 0°, 7°, 12° and 18° Lordosis\*

\* Some sizes and angles available by request

Amplify Surgical is a privately-held spinal device company that aims to revolutionize **minimally-invasive spine surgery** by transforming ordinary procedures with groundbreaking surgical solutions.

For more information or to place an order:

**765-AMPLIFY (765-267-5439) | [info@amplifysurgical.com](mailto:info@amplifysurgical.com) | [amplifysurgical.com](http://amplifysurgical.com)**



Please refer to package insert for current warnings, precautions, and instructions for use.