



dualX transforms the fusion environment by delivering a powerful implant through a minimally-invasive approach

Please send your completed scrub sheets to:  
[customerservice@amplifysurgical.com](mailto:customerservice@amplifysurgical.com)  
 765-AMPLIFY (765-267-5439)

Top View



Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
<b>Lateral Lumbar Interbody Fusion (LLIF)</b>								
1283-1005-002	dualX LLIF 7-9mm x 45mm	0°	13	22	7	9	54	45
1283-1005-003	dualX LLIF 7-9mm x 50mm	0°	13	22	7	9	59	50
1283-1005-004	dualX LLIF 7-9mm x 55mm	0°	13	22	7	9	64	55
1283-1007-001	dualX LLIF 8-11mm x 40mm	7°	13	22	8	11	49	40
1283-1007-002	dualX LLIF 8-11mm x 45mm	7°	13	22	8	11	54	45
1283-1007-003	dualX LLIF 8-11mm x 50mm	7°	13	22	8	11	59	50
1283-1007-004	dualX LLIF 8-11mm x 55mm	7°	13	22	8	11	64	55
1283-1007-012	dualX LLIF 10-13mm x 45mm	7°	13	22	10	13	54	45
1283-1007-013	dualX LLIF 10-13mm x 50mm	7°	13	22	10	13	59	50
1283-1007-014	dualX LLIF 10-13mm x 55mm	7°	13	22	10	13	64	55
1283-1007-022	dualX LLIF 12-15mm x 45mm	7°	13	22	12	15	54	45
1283-1007-023	dualX LLIF 12-15mm x 50mm	7°	13	22	12	15	59	50
1283-1007-024	dualX LLIF 12-15mm x 55mm	7°	13	22	12	15	64	55
1283-1000-006	dualX LLIF 10-13mm x 40mm	12°	13	22	10	13	49	40
1283-1000-007	dualX LLIF 10-13mm x 45mm	12°	13	22	10	13	54	45
1283-1000-008	dualX LLIF 10-13mm x 50mm	12°	13	22	10	13	59	50
1283-1000-009	dualX LLIF 10-13mm x 55mm	12°	13	22	10	13	64	55
1283-0009-001	dualX LLIF Lock Screw 2 - 40mm					N/A		
1283-0009-002	dualX LLIF Lock Screw 2 - 45mm					N/A		
1283-0009-003	dualX LLIF Lock Screw 2 - 50mm					N/A		
1283-0009-004	dualX LLIF Lock Screw 2 - 55mm					N/A		

Note: dualX expansion in width and height is definitive and not variable. For example, the dualX LLIF 0 DEG x 7-9mm x 45mm expands exactly from 13mm to 22mm (9mm expansion) in width and 7mm to 9mm (2 mm) in height.

Please refer to package insert for current warnings, precautions and instructions for use.  
 Copyright 2025. All Rights Reserved. dualX is a registered trademark of Amplify Surgical, Inc.  
 For more information call 765-AMPLIFY, email [info@amplifysurgical.com](mailto:info@amplifysurgical.com) or visit [amplifysurgical.com](http://amplifysurgical.com)  
 ML.10014.G

