

Top View




Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal Lumbar Interbody Fusion (TLIF)								
1253-1006-001	dualX TLIF 7-9mm	0°	12	21	7	9	36	30
1253-1006-004	dualX TLIF 8-11mm	8°	12	21	8	11	36	30
1253-1006-005	dualX TLIF 9-12mm	8°	12	21	9	12	36	30
1253-1006-006	dualX TLIF 10-13mm	8°	12	21	10	13	36	30
1253-1006-007	dualX TLIF 11-14mm	8°	12	21	11	14	36	30
1253-1006-008	dualX TLIF 12-15mm	8°	12	21	12	15	36	30
1253-0033	dualX TLIF Lock Screw 2		N/A					

Top View



Side View



Part Number	Description	Lordosis	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal/Posterior Lumbar Interbody Fusion (T/PLIF)								
1282-1000-001	dualX T/PLIF 7-9mm	0°	12	21	7	9	36	22
1282-1001-014	dualX T/PLIF 8-10mm	8°	12	21	8	10	36	22
1282-1000-004	dualX T/PLIF 8-11mm	8°	12	21	8	11	36	22
1282-1000-005	dualX T/PLIF 9-12mm	8°	12	21	9	12	36	22
1282-1000-006	dualX T/PLIF 10-13mm	8°	12	21	10	13	36	22
1282-1000-007	dualX T/PLIF 11-14mm	8°	12	21	11	14	36	22
1282-1001-016	dualX T/PLIF 10-12mm	12°	12	21	10	12	36	22
1282-1001-017	dualX T/PLIF 10-13mm	12°	12	21	10	13	36	22
1282-1001-018	dualX T/PLIF 11-13mm	15°	12	21	11	13	36	22
1282-1000-019	dualX T/PLIF 11-14mm	15°	12	21	11	14	36	22
1282-0005	dualX T/PLIF Lock Screw 2		N/A					